

**MP1651-2-3-4 PASTES FOR IMPREX
RED, YELLOW, BLACK, ORANGE**

Product : Paste, free flowing or slightly thixotropic.

Introduction: The PASTA X IMPREX have not only a maximum transparency but also a beneficial influence on the resistance to weathering as a result of their high UV absorption.

Processing : Colour preparations are best incorporated into the 2873 by stirring vigorously. Should speckles develop due to a rapid change of concentration, a two-step process may be preferred: 1. prepare a concentrated product (2873 10% + 4 – 6% MP151-2-3-4), 2. incorporate in the remaining 2873.

Available colours: MP1651 pasta x imprex Red, MP1652 pasta x imprex yellow, MP1653 pasta x imprex black, MP1654 pasta x imprex Orange.

Storage: If kept in unopened original containers and at temperatures below 30°C, the product can be stored for 8 months. Once opened, a container's content should be used up within a short time. When stored at temperatures above 30°C, evaporation of the solvent and/or chemical reaction of the binder may produce a surface skin in PASTA X IMPREX. After completely removing the skin, the remaining content is fit for use. Under no circumstances should the skin be dispersed by stirrer or dissolver.

General notes : Always mix the products before use